




COOLPALL™ VERTOS ADVANCE



OUR NEW BULK SHIPPING SYSTEM FEATURING PCM COOLANTS

CoolPall™ Vertos Advance is the NEW high performance bulk shipping system from Pelican BioThermal™. This bulk shipper features PCM coolants and has been designed to be flat packed, modular, and easily assembled. The CoolPall™ Vertos Advance maintains 2°C to 8°C for 120 hours of temperature controlled performance.

WHY CHOOSE COOLPALL™ VERTOS ADVANCE?

Ease & Speed of Use

- Superior performance foam and PCM technology reduces excursions and allows for hibernation to be achieved during transit, which pauses phase change within the shipper and can yield longer temperature compliant shipments
- Ergonomic design utilizes innovative interlocking components to ensure quick and easy assembly
- Modular components reduces inventory and allows for simple ordering of spare parts
- Reduced flat-packed height minimizes in-bound distribution costs and frees up valuable warehouse storage space
- Accepts full pallets of product for operational speed and simplicity

Robust performance

- Provides strict temperature control of product within the demanding environment of global distribution
- Utilizes high density moulded materials for thermal and mechanical robustness while remaining lightweight
- Best in class payload to external volume ratio, lowering transport costs

COOLPALL™ VERTOS ADVANCE

DURATION

120+ Hours

TEMPERATURE

2°C to 8°C



VOLUME

ISO Pallet size,
volume up to 1680L

COOLPALL™ VERTOS ADVANCE PRODUCT SPECIFICATION

Product	Exterior Dimensions - Inches	Interior Dimensions - Inches	Payload Volume	Tare Weight - lbs
ISOPAD	53 x 62 x 63	44 x 51 x 45	1680L	518

Peli BioThermal Profile, CoolPall™ Vertos Advance +2°C to +8°C



Interlocking high density construction
Increases the possibility of re-use



Inspection hatch
Minimizing the potential for excursions during customs inspection



FOR MORE INFORMATION ABOUT PELICAN BIOTHERMAL PRODUCTS, PLEASE CONTACT:

INFO@PELICANBIOTHERMAL.COM T: +1 763 412 4800 PELICANBIOTHERMAL.COM

ISO17025, ISO9001-2008 Accredited

Golden Hour™, TIC™, Crêdo™, Crêdo Cube™, Crêdo ProEnvision™, Crêdo Xtreme™, Crêdo DURACUBE™, Crêdo ProMed™, Sherpa™ Systems, Chronos™ Express, Chronos™ Advance and CoolPall™ Vertos are registered trademarks of Pelican Products, Inc., its subsidiaries and/or affiliates.



©2018 PELICAN BioThermal all rights reserved.