

Credo Duracube™ Product Line



DURACUBE LT

SERIES 4 / 22 / 20M

16 Liters - 96 Hours

DURACUBE HD

SERIES 4

2 Liters - 36 Hours

16 Liters - 96 Hours

28 Liters - 96 Hours

56 Liters - 168 Hours

96 Liters - 96 Hours

SERIES 22

16 Liters - 96 Hours

28 Liters - 96 Hours

56 Liters - 168 Hours

96 Liters - 96 Hours

SERIES 20M

16 Liters - 96 Hours

28 Liters - 96 Hours

56 Liters - 168 Hours

96 Liters - 96 Hours

SERIES 50M

12 Liters - 96 Hours

56 Liters - 96 Hours

96 Liters - 72 Hours

EXTREMELY DURABLE HARD SHELL THERMALLY PROTECTED SHIPPERS

The Crēdo Duracube™ product line offers a selection of superior thermal protection in a comprehensive range of hard shell outer containers (between 2 liters and 96 liters) designed to strike an optimal balance between case weight, durability and duration of transit capabilities while providing a 10-15% increase in thermal performance over corrugated material.



The Crēdo DURACUBE™ LT (high density polyethylene) is qualified to hold room-temperature (CRT), refrigerated and frozen (25°C to -40°C) medical materials such as platelets and bio-pharma products at a safe temperature for up to four days (96 hours).

- Longer life span, lower cost of ownership and less maintenance expense.
- Reusable and lightweight durability - only 2-3 lbs. more tare weight than plastic corrugate.
- Ideal for use in reducing shipping costs and eliminating payload excursions when unexpected delays may be eminent.
- Modular, beveled TIC™ walls with integrated phase-change material break down easily for simple storage and preconditioning. VIP's are individually replaceable.
- Single packout assembly for all temperature ranges provides ease of use and simplified training requirements - less distribution costs and lower cost per per use.



Shipper Name	Payload Interior Dimensions (inches)	Outer Dimensions (inches)	Tare Weight (lbs.)	ISTA 7D Summer / Winter Performance	Volume (liters)
Series 4 1696 LT	10 x 10 x 10	17.5 x 17.5 x 16.25	34.8	93.42 / 138.67	16
Series 22 1696 LT	10 x 10 x 10	17.5 x 17.5 x 16.25	34.8	126.5 / 68.25	16
Series 20M 1696 LT	10 x 10 x 10	17.5 x 17.5 x 16.25	40.6	97.75 / X	16
Series 4 2896 LT	12 x 12 x 12	19.6 x 19.6 x 17.8	48.0	98.8 / 135.75	28

SERIES 4 (2-8°C)

GH0236
12168
1696
2896
56168
9696

SERIES 22 (15-25°C)

1696
2896
56168
9696



SERIES 20M (-18°C and colder)

1696
2896
56168
9696

SERIES 50M (-40°C and colder)

1296
5696
9672



Shipper Name	Payload Interior Dimensions (inches)	Outer Dimensions (inches)	Tare Weight (lbs.)	ISTA 7D Summer Performance	ISTA 7D WINTER Performance	Volume (liters)
S4 GH0236	6 x 5 x 4.25	11.8 x 9.8 x 7.7	9.7	41.17	75.25	2
S4 12168 HD	8.9 x 8.9 x 8.9	20.5 x 20.5 x 19.25	55.5	169.92	170+	12
S4 1696 HD	8.90 x 8.9 x 8.9	20.5 x 20.5 x 19.25	45.20	102.00	85.83	16
S4 2896 HD	12 x 12 x 12	21 x 21 x 19.25	53.8	106.33	187.83	28
S4 56168 HD	15 x 15 x 15	26.5 x 26.5 x 25.25	105.0	192.25	235+	56
S4 9696 HD	18 x 18 x 18	30 x 25.75 x 28	105.2	128.83	144+	96
S22 1696 HD	10 x 10 x 10	20.5x 20.5 x 19.25	45.0	102.00	85.83	16
S22 2896 HD	12 x 12 x 12	21 x 21 x 19.25	X	141.42	145.58	28
S22 56168 HD	15 x 15 x 15	26.5 x 26.5 x 25.25	105.2	192.00	147.33	56
S22 9696 HD	18 x 18 x 18	27 x 27 x 26	105.2	142.83	114.0	96
S20M 1696 HD	10 x 10 x 10	21 x 21 x 19.25	51.0	118.25	X	16
S20M 2896 HD	12 x 12 x 12	21 x 21 x 19.25	62.2	120.92	X	28
S20M 56168 HD	15 x 15 x 15	26.5 x 26.5 x 25.25	116.8	123.33	X	56
S20M 9696 HD	18 x 18 x 18	26.5 x 26.5 x 26	123.0	139.75	X	96
S50M 1296 HD	8.9 x 8.9 x 8.9	18 x 16.5 x 17	38.0	71.33	X	12
S50M 5696 HD	15 x 15 x 15	26.5x 26.5 x 25.25	121.0	98.08	X	56

PELICAN BIOTHERMAL, LLC
Sales & Operations Headquarters
3020 Niagara Lane North
Plymouth, MN 55447 USA
1.877.537.9800
www.pelicanbiothermal.com
B7_V6_0915

All trademarks are registered and/or unregistered trademarks
of Pelican BioThermal LLC, its parent, subsidiaries and/or
affiliates.

ISO 9001-2008 Certification
ISO14001-2004 Certification
OHASA18001-2007 Certification
ISO17025 Certification

