




Crēdo™

ON DEMAND



THE NEW HIGH PERFORMANCE, COST EFFECTIVE RENTAL PROGRAMME FOR TEMPERATURE CONTROLLED PALLET SHIPPERS



WHY CHOOSE CRĒDO™ ON DEMAND?

- Durable passive temperature-controlled system that is less costly than active systems, easy implementation
- Global availability, convenient pick-up and drop-off locations
- Pre-conditioning available for quick 'load and go'
- Flexible daily rentals available (minimum 6 days)
- Two options available, Crēdo Xtreme™ for volumes upto 807L, or Crēdo™ Cargo for volumes up to 1686L

FOR MORE INFORMATION ABOUT PELICAN BIOTHERMAL PRODUCTS, PLEASE CONTACT:
INFO@PELICANBIOTHERMAL.COM T: +1 763 412 4800 PELICANBIOTHERMAL.COM

ISO17025, ISO9001-2008 Accredited

Golden Hour™, TIC™, Crēdo™, Crēdo Cube™, Crēdo ProEnvision™, Crēdo Xtreme™, Crēdo Duracube™, Crēdo ProMed™, Sherpa Systems™, CoolGuard™ and CoolPall™ are registered trademarks of Pelican BioThermal LLC., its subsidiaries and/or affiliates.

©2020 PELICAN BioThermal all rights reserved.

TEMPERATURE RANGE

2°C to 8°C Chilled
15°C to 25°C Ambient
-20°C Frozen

DURATION

120 - 144+ hours
(ISTA 7D)

CRĒDO XTREME™

OUTER DIMENSIONS
1226 x 1022 x 1340mm

PAYLOAD DIMENSIONS
1027 x 825 x 953mm

VOLUME
807 LITERS

CRĒDO™ CARGO

OUTER DIMENSIONS
1560 x 1525 x 1600mm

PAYLOAD DIMENSIONS
1295 x 1085 x 1200mm

VOLUME
1686 LITERS