



# Crēdo DURACUBE™



## DURABLE HARD SHELL THERMALLY PROTECTED SHIPPERS

The Crēdo DURACUBE™ product line offers a selection of superior thermal protection in a comprehensive range of lightweight and heavy duty hard shell outer containers (between 2 liters and 96 liters) designed to strike an optimal balance between case weight, durability and duration of transit capabilities while providing a 10-15% increase in thermal performance over outer corrugated material.



## CRĒDO DURACUBE™ LT

The Crēdo DURACUBE™ LT (high density polyethylene) is qualified to hold room-temperature (CRT), refrigerated and frozen (25°C to -40°C) medical materials such as platelets and bio-pharma products at a safe temperature for up to four days (96 hours).

- Longer life span, lower cost of ownership and less maintenance expense
- Reusable and lightweight durability - only 2-3 lbs. more tare weight than plastic corrugate
- Ideal for use in reducing shipping costs and eliminating payload excursions when unexpected delays may be imminent
- Modular, beveled TIC™ walls with integrated phase-change material break down easily for simple storage and conditioning. VIP's are individually replaceable
- Single packout assembly for all temperature ranges provides ease of use and simplified training requirements - less distribution costs and lower cost per per use



### CRĒDO DURACUBE™ LT

**DURATION**  
120+ Hours

**SERIES 4** (2°C to 8°C)

**VOLUME**  
16L, 28L

**SERIES 22** (15°C to 25°C)

**VOLUME**  
16L

**SERIES 20M** (-25°C to -15°C)

**VOLUME**  
16L

### CRĒDO DURACUBE™ HD

**DURATION**  
120+ Hours

**SERIES 4** (2°C to 8°C)

**VOLUME**  
2L, 16L, 28L, 56L, 96L

**SERIES 22** (15°C to 25°C)

**VOLUME**  
16L, 28L, 56L, 96L

**SERIES 20M**  
(-18°C and colder)

**VOLUME**  
16L, 28L, 56L, 96L

**SERIES 50M**  
(-40°C and colder)

**VOLUME**  
12L, 56L, 96L

## CRĒDO DURACUBE™ HD

The CrĒdo DURACUBE™ HD (HPX™ resin exterior) is highly durable and able to secure thermal protection and endurance in rigorous conditions. Coupled with our patented TIC™ system with encapsulated PCM and surrounded by vacuum insulated panels, the CrĒdo DURACUBE™ HD product line protects critical payloads for up to seven days at a safe temperature (168 hours).

- 10-15% increase in thermal performance over corrugated material
- Designed and tested to accommodate transport of fragile, expensive, and larger pharmaceutical shipments
- HD systems offer safe temperature ranges from 25°C to -40°C and colder
- Ideal for mobile blood bank use, overseas shipping, or situations requiring rugged transport of delicate medical materials



## CRĒDO DURACUBE™ LT

Type	Liters	Interior Dimensions - Inches	Exterior Dimensions - Inches	Tare Weight -lbs
Series 4 1696 LT	16L	10 x 10 x 10	18 x 18 x 16	33
Series 22 1696 LT	16L	10 x 10 x 10	18 x 18 x 16	33
Series 20M 1696 LT	16L	10 x 10 x 10	18 x 18 x 16	39
Series 4 2896 LT	28L	12 x 12 x 12	20 x 20 x 18	50

## CRĒDO DURACUBE™ HD

Type	Liters	Interior Dimensions - Inches	Exterior Dimensions - Inches	Tare Weight -lbs
<b>SERIES 4</b>				
S4 GH0236	2L	6 x 5 x 4	12 x 10 x 8	10
S4 12168 HD	12L	9 x 9 x 9	21 x 21 x 19	57
S4 1696 HD	16L	10 x 10 x 10	21 x 21 x 19	47
S4 2896 HD	28L	12 x 12 x 12	21 x 21 x 19	57
S4 56168 HD	56L	15 x 15 x 15	27 x 27 x 25	105
S4 9696 HD	96L	18 x 18 x 18	27 x 27 x 25	102
<b>SERIES 22</b>				
S22 1696 HD	16L	10 x 10 x 10	21 x 21 x 19	47
S22 56168 HD	56L	15 x 15 x 15	27 x 27 x 25	105
S22 9696 HD	96L	18 x 18 x 18	27 x 27 x 25	102
<b>SERIES 20M</b>				
S20M 1696 HD	16L	10 x 10 x 10	21 x 21 x 19	53
S20M 2896 HD	28L	12 x 12 x 12	21 x 21 x 19	66
S20M 56168 HD	56L	15 x 15 x 15	27 x 27 x 25	117
S20M 9696 HD	96L	18 x 18 x 18	27 x 27 x 25	119
<b>SERIES 50M</b>				
S50M 5696 HD	56L	15 x 15 x 15	27 x 27 x 25	119
S50M 9672 HD	96L	18 x 18 x 18	27 x 27 x 25	124

FOR MORE INFORMATION ABOUT PELICAN BIOTHERMAL PRODUCTS, PLEASE CONTACT:

INFO@PELICANBIOTHERMAL.COM T: +1 763 412 4800 PELICANBIOTHERMAL.COM

ISO17025, ISO9001-2008 Accredited

Golden Hour™, TIC™, CrĒdo™, CrĒdo Cube™, CrĒdo ProEnvision™, CrĒdo Xtreme™, CrĒdo DURACUBE™, HPX™, CrĒdo ProMed™, Sherpa™ Systems, Chronos™ Express, Chronos™ Advance and CoolPall™ Vertos are registered trademarks of Pelican Products, Inc., its subsidiaries and/or affiliates.



©2017 PELICAN BioThermal all rights reserved.