



CRĒDO CUBE™ SERIES 20M & 50M

Superior thermal protection in a convenient mid-size container.

The Crēdo Cube™ is a passive and reusable shipper qualified to hold chilled medical materials at a safe temperature for up to five days. This system is ideal for reducing shipping costs or for circumstances in which unexpected delays may occur.

Modular TIC™ (Thermal Isolation Chamber) walls with integrated PCM (Phase-Change Material) and VIPs (Vacuum Insulated Panels) makes conditioning in a freezer for year-round shipping easy and allows effortless packing and storage.

WHY CHOOSE CRĒDO™ SYSTEMS?

- The most trusted and well recognized brand of high-performance reusable shippers with the largest number of sizes and temperature ranges available
- Unique TIC™ walls contain PCM and completely enclose the temperature-sensitive payload ensuring consistently stable performance
- Crēdo™ shippers are highly durable, lightweight, and provide an optimum volume ratio of interior payload space for your products - allowing you to ship more product in fewer containers, reducing distribution costs



ENHANCE YOUR CRĒDO™ SYSTEMS

- **Crēdo ProEnvision™**: Easy-to-use, secure, web-based asset tracking software and global service centers available to assist customers and maximize their return on investment

CRĒDO CUBE™

DURATION

96 - 168 Hours

SERIES 20M

(-25°C to -15°C)

VOLUME

3L, 4L, 8L, 10L, 12L, 16L, 28L, 56L

SERIES 50M

(-40°C and Cooler)

VOLUME

4L, 12L, 16L, 28L

CRĒDO™ SERIES 20M -25°C TO -15°C

Payload Size	Exterior Dimensions - Inches	Interior Dimensions - Inches	Tare Weight - lbs
Single Wall VIP, White Plastic Corrugate			
3L	13 x 11 x 11	5 x 5 x 5	16
4L	6 x 6 x 6	6 x 6 x 6	14
8L	18 x 11 x 12	11 x 7 x 7	20
10L	15 x 14 x 14	9 x 9 x 9	30
12L	15 x 14 x 14	9 x 9 x 9	27
16L	17 x 15 x 15	10 x 10 x 10	31
28L	18 x 17 x 17	12 x 12 x 12	45
Single Wall VIP, Kraft Box			
56L	21 x 20 x 20	15 x 15 x 15	66
Nested VIP, White Plastic Corrugate			
4L	15 x 14 x 14	6 x 6 x 6	22
12L	18 x 17 x 17	9 x 9 x 9	41
Nested VIP, Kraft Box			
28L	21 x 20 x 20	12 x 12 x 12	65
Double Wall VIP, White Plastic Corrugate			
16L	18 x 17 x 18	10 x 10 x 10	43

CRĒDO™ SERIES 50M (-40°C AND COOLER)

Payload Size	Exterior Dimensions - Inches	Interior Dimensions - Inches	Tare Weight - lbs
Nested VIP, White Plastic Corrugate			
4L	15 x 14 x 14	6 x 6 x 6	22
12L	18 x 17 x 17	9 x 9 x 9	42
Nested VIP, Kraft Box			
28L	21 x 20 x 20	12 x 12 x 12	67
Double Wall VIP, White Plastic Corrugate			
16L	18 x 17 x 18	10 x 10 x 10	44



Superior Performance
ensuring complete
product integrity

Small and Light
significantly reducing
distribution costs

**Crêdo ProEnvision™ and
Service Centers**
support reclamation
and re-use

FOR MORE INFORMATION ABOUT PELICAN BIOTHERMAL PRODUCTS, PLEASE CONTACT:

INFO@PELICANBIOTHERMAL.COM T: +1 763 412 4800 PELICANBIOTHERMAL.COM

ISO17025, ISO9001-2008 Accredited

Golden Hour™, TIC™, Crêdo™, Crêdo Cube™, Crêdo ProEnvision™, Crêdo Xtreme™, Crêdo DURACUBE™, Crêdo ProMed™, Sherpa™ Systems, Chronos™ Express, Chronos™ Advance and CoolPall™ Vertos are registered trademarks of Pelican Products, Inc., its subsidiaries and/or affiliates.



©2017 PELICAN BioThermal all rights reserved.